

Pushing Standards of Infusion Therapy Even Higher

ambIT® zeros in on the most essential features to ensure patient's needs are met in a **simple, safe, and adaptable** manner.

ambIT® PCA:

Patient controlled bolus pump that maximizes functionality while eliminating cost and waste of disposable devices.

ambIT® PreSet:

Faster and easier programming with four pain management protocols to choose from and a fifth option to customize your own.



Better by design, the ambIT® pump is a safe, simple, multipurpose fluid delivery system

Flexibility and Convenience:

- Small and lightweight design
- Easily concealed in custom-built carry pouch with view-window
- Easy-press bolus button
- Titratable infusion rates
- History reports
- No maintenance

Effective Safety Standards:

- Programming lockout
- Bolus lockout
- Anti-free flow
- Program review
- Device internal self-check
- Latex-free and DEHP-free fluid path

504 West 8360 South
Sandy, Utah 84070 | USA
801-352-1888
877-352-1888
Fax: 801-352-1818
service@ambitpump.com

ambIT® is a division of
Summit Medical Products, Inc.



Volumetric Accuracy	+/- 6%
Basal Flow Rates	0.1 to 20 ml/hr (increments vary by model)
Stroke Volume	50 microliters
Bolus Delivery Rate	100 ml/hr
Bolus Dose Volume (ml)	0.1 to 20 ml
Bolus Lockout Times (hr:min)	00:01 to 12:00 (increments vary by model)
Volumes to be Infused (ml)	25 to 1000 ml
Pump Mechanism	Microprocessor Controlled Rotary-Peristaltic
Dimensions	2.16 in. x 1.4 in. x 6.875 in. (55 mm x 36 mm x 175 mm)
Weight	6.4 ounces (181.4 grams) With Batteries
Power Supply	Two AA Batteries
Alarms & Alerts	<ul style="list-style-type: none"> • Detached Cassette • Malfunction • Run Indicator Light • Bolus Infusing • Occlusion Downstream • Low Battery • Dead Battery • Boluses Requested • Infusion Complete • Pause Indicator
Reports	<ul style="list-style-type: none"> • Volume Infused • Boluses Requested • Boluses Delivered • Time Elapsed
Delivery Profiles	Basal Rate + Bolus Basal Rate Only Bolus Only

The ambIT® Pump follows internationally recognized volumetric accuracy standards (IEC 60601-2-24), unlike elastomeric and mechanical pumps. Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

504 West 8360 South
Sandy, Utah 84070 | USA
801-352-1888
877-352-1888

Fax: 801-352-1818
service@ambitpump.com

ambIT® is a division of
Summit Medical Products, Inc. 