

Designed Specifically for Higher Infusion Rates

ambIT® Continuous can run at higher infusion rates without the option of a bolus - essential for applications in chemotherapy, antibiotic treatments, and hydration services.

Durable and Reusable:

ambIT® Continuous is designed to eliminate mechanical maintenance and calibration expenses.

Patients can easily carry the ambIT® Continuous while receiving accurate infusion therapy.



Effective Safety Standards:

- Programming lockout
- Anti-free flow
- Program review
- Device internal self-check
- Latex-free and DEHP-free fluid path

Flexibility and Convenience:

- Small and lightweight design
- Easily concealed in custom-built carry pouch with view-window
- Titratable infusion rates
- Infusion history reports
- No maintenance

Only Pertinent Alarms:

- Downstream Occlusion Alarm
- Check Cassette Alarm
- Product Malfunction

504 West 8360 South
Sandy, Utah 84070 | USA
801-352-1888
877-352-1888
Fax: 801-352-1818
service@ambitpump.com

ambIT® is a division of
Summit Medical Products, Inc.



ambIT® CONTINUOUS

Volumetric Accuracy	+/- 6%
Flow Rates	0.1 to 125 ml/hr
Stroke Volume	50 microliters
Volumes to be Infused (ml)	25 to 999 ml in 1 ml increments
Pump Mechanism	Microprocessor Controlled Rotary-Peristaltic
Dimensions	2.16 in. x 1.4 in. x 6.875 in. (55 mm x 36 mm x 175 mm)
Weight	6.4 ounces (181.4 grams) With Batteries
Power Supply	Two AA Batteries
Alarms & Alerts	<ul style="list-style-type: none">• Detached Cassette• Malfunction• Run Indicator Light• Occlusion Downstream• Low Battery• Dead Battery• Infusion Complete• Pause Indicator
Reports	Volume Infused Time Elapsed
Delivery Profiles	Continuous Rate

The ambIT® Continuous Pump follows internationally recognized volumetric accuracy standards (IEC 60601-2-24), unlike elastomeric and mechanical pumps. Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

504 West 8360 South
Sandy, Utah 84070 | USA
801-352-1888
877-352-1888

Fax: 801-352-1818
service@ambitpump.com

ambIT® is a division of
Summit Medical Products, Inc. 

40069.00 1/16

www.ambitpump.com